

NEW

unoQUIP

CATHETERISATION TRAY WITH URINE METER

Complete Foley Catheter Kit



1

Complete Catheterisation Tray

Supports aseptic
catheterisation and an
efficient clinical workflow

2

Preconnected 100% Silicone Foley Catheter

Designed for smooth
catheterisation and reliable
drainage

3

Integrated Q500 Plus Diuresis Measuring System

Safe and accurate urine
measurement in critically ill
patients

Complete Foley Catheter Kit

The Foley Catheter Kit combines all essential components to support aseptic catheterisation, reliable urinary drainage and safe urine measurement.

1 Complete Catheterisation Tray



Catheterisation procedure supplies:

- 1 x Waste Bag (orange)
- 1 x Apron
- 2 x Absorbent Underpads
- 1 x Fenestrated Drape
- 2 x Pairs of Nitrile Gloves
- 1 x 10% Aqueous Glycerine Solution
- 2 x Cleansing Solutions
- 1 x Water Based Lubricant
- 4 x Gauze Pads 10 x 10 cm
- 1 x Empty Syringe
- Peel-off Labels for full traceability and documentation of catheterisation details

2 Preconnected 100% Silicone Foley Catheter



- Soft, flexible silicone with wide drainage lumen to reduce blockage risk
- Long-term usage for up to 90 days
- Pre-attached to urine meter with tamper-evident seal to help prevent disconnection

3 Integrated Q500 Plus Diuresis Measuring System



- Top-mounted system with two Non-Return Valves to reduce infection risk
- Double-Lumen Free-Flow System for continuous urine flow and complete drainage
- Designed for efficient clinical handling and patient transport

ORDERING INFORMATION

NHS Ref.	UNOQUIP Ref.	Catheter size	Pack QTY.
FUU85099	T01-41151210-01EN	CH 12	1/5
FUU85098	T01-41151410-01EN	CH 14	1/5
FUU85097	T01-41151610-01EN	CH 16	1/5
FUU85096	T01-41151810-01EN	CH 18	1/5

All configurations include:

Q500 Plus Diuresis Measuring System with 1.5 m drainage tube, 2 L bag, hook and hanging straps; 100% Silicone Foley Catheter with 10 ml balloon; catheterisation procedure supplies.

CONTACT US



Scan to explore
product details
on our website

📍 **UNOQUIP UK**
27 Old Gloucester Street
London WC1N 3AX
United Kingdom

✉ uksupport@unoquip.com

📍 **UNOQUIP GmbH**
Schützengraben 20
8200 Schaffhausen
Switzerland

🌐 www.unoquip.com

